
Temperature Controllers

CTC100 — Cryogenic temperature controller



- 4 temperature sensor inputs
- 2 powered & 4 analog voltage outputs
- Up to 6 feedback control loops
- 4 analog & 8 digital I/O channels
- Graphical touchscreen display
- Data logging on removable flash media
- User programs (macros)
- USB, Ethernet, RS-232 & GPIB (opt.)

• CTC100 ... \$2600 (U.S. list)

CTC100 Temperature Controller

Introducing the new CTC100 Cryogenic Temperature Controller — a high performance instrument that can monitor and control temperatures with millikelvin resolution.


The CTC100 Cryogenic Temperature Controller is configured to suit a wide range of research and industrial applications. The system consists of four sensor inputs, two powered and four analog voltage outputs, and up to six feedback control loops. Four general purpose analog and eight digital I/Os are available, along with autotuning functions for setting PID parameters automatically.

Sensor Inputs

The CTC100 offers four temperature inputs that can read RTDs, thermistors, and diodes. Each temperature input channel has its own 24-bit ADC with eleven input ranges, and is equipped with its own independent excitation current source.

Standard calibration curves for a variety of sensors are included, and custom calibration curves of up to 200 points each can be entered. Each sensor input has high and low level or rate-of-change alarms. Sensor inputs can be lowpass-filtered to reduce noise, and/or differenced with another channel.





Channel setup menu

Powered and Unpowered Outputs

The CTC100 has two heater outputs that can each deliver up to 100 W of power to a 25 Ω heater. In addition, four analog voltage I/O channels can be used to drive heaters with the help of an external amplifier.

PID Feedback Control

With up to six feedback control loops available, the CTC100 can provide precise temperature control of each of its heater outputs by continually adjusting the heater power. Any of the CTC100's channels can be selected as the input for each feedback loop. Feedback time constants can be adjusted between 200 ms and 10 hrs.

Up to ten sets of PID parameters can be stored for each channel. Setpoints can be ramped at a fixed rate (or with a user program), set from an analog input.

Analog and Digital I/Os

The CTC100 comes with four general-purpose ± 10 V voltage I/O channels read by a 24-bit ADC. It also has eight digital I/O channels that can interact with user programs.

Four 5 A relays can be used for process control. Three virtual channels, not connected to any physical input, allow calculated values (such as the difference between two channels, or a value calculated by a user program) to be displayed, graphed, and logged.

Numeric and Graphical Display

The CTC100's color LCD display can show any combination of temperature measurements and heater outputs on graphs or numeric displays. Up to eight channels can be plotted either on a single graph with a common Y axis, or on separate graphs with independent Y axes. Touchscreen operation makes the instrument versatile and easy to use.


1 Single	Multiple	Custom	Ponytail	Numerio
ln 1			50.1	178°
ln 3			25.7	773°
ln 4			25.8	311°

Numeric display

Data Logging

Up to 4096 readings per channel can be logged to the CTC100's internal memory. For longer-term storage, data can be logged to standard USB memory sticks or hard drives.

Data logged to USB devices can be transferred to a computer by plugging the USB device into a PC. Windows applications are included to graph CTC100 log files and to convert them to various ASCII text formats.



Graphical display



Computer Communications

Each of the CTC100's front-panel controls has a corresponding text command that can be sent over USB, Ethernet, and either RS-232 or an optional GPIB interface.

When the USB interface is used, the CTC100 appears on the computer as a standard COM port and can be controlled by any software that is compatible with an RS-232 port.



Select menu

User Programs

User programs (macros), consisting of one or more remote commands, can be uploaded to the CTC100. This is done by either sending them through one of the communications ports or by saving them as text files on a USB memory device and then plugging the device into the CTC100. Program macros can also be entered and edited from the front panel.



CTC100 rear panel

Specifications

Temperature Controller

Min. sampling rate 1 Hz Max. sampling rate 50 or 60 Hz

(depending on AC line frequency)

Data logging rate 10 samples/second/channel to

1 sample/hour/channel

(can be set independently for each channel or globally for all channels) 0.001 °C, °F, K, V, A, W, etc.

Display resolution $0.001 \,^{\circ}\text{C}$, $^{\circ}\text{F}$, K, V, A, W, etc. if -1000 < displayed value < 1000;

6 significant figures otherwise

PID autotuning Single step response or relay tuning with conservative, moderate, and

aggressive response targets

Display 320×240 pixel color touchscreen; numeric and graphical data displays. Alarms Upper and lower temperature limits

or rate-of-change limits can be

set on each channel.

Computer interface USB, Ethernet, and RS-232;

GPIB (IEEE488.2) optional 10 A, 88 to 132 VAC

Power 10 A, 88 to 132 VAC or 176 to 264 VAC, 47 to 63 Hz

or DC

Dimensions $8.5" \times 5" \times 16"$ (WHL)

Weight 13 lbs.

Warranty One years parts and labor on defects

in material and workmanship

Analog I/O

Inputs/outputs 4 voltage I/O channels, independently

configurable as inputs or outputs

Connector 4 BNC jacks Range $\pm 10 \text{ V}$

Resolution 24-bit input, 16-bit output ADC noise 30 µVrms (at 10 samples/s)

Digital I/O

Inputs/outputs 8 optoisolated TTL lines,

configurable as either 8 inputs or

8 outputs

Connector One DB-25F

Relays

Outputs 4 independent SPDT relays

100 W DC outputs

Output Two unipolar DC current sources

Connector #6 screw terminals

Range 50 V 2 A, 50 V 0.2 A, 50 V 0.02 A,

20 V 2 A, 20 V 0.2 A, 20 V 0.02 A



CTC100 Specifications

Output resolution 16 bit

Accuracy $\pm 1 \text{ mA } (2 \text{ A range})$

 $\pm 0.02 \,\text{mA} \, (0.2 \,\text{A range}) \\ \pm 0.002 \,\text{mA} \, (0.02 \,\text{A range})$

Noise (rms) $(25 \Omega \text{ load, DC to } 10 \text{ Hz})$

5 μA (2 A range) 0.5 μA (0.2 A range) 0.05 μA (0.02 A range)

Diodes, Thermistors and RTD inputs

Inputs Four inputs for 2-wire or 4-wire thermistor, diode, or RTD

Socket Two DB9 (female)

	()					
	Input Range	Excitation Current	Initial Accuracy	Temp. Drift (typ.) (at midrange)	Noise (rms) (at midrange)	
Diodes	0 to 2.5 V	10 μΑ	$10 \mu\text{V} + 0.01 \%$ of rdg	±5 ppm/°C	3 μV	
RTDs	0 to 10Ω	3 mA	±0.005 Ω	±0.0001 Ω/°C	0.0001Ω	
	0 to 30Ω	3 mA	$\pm 0.005\Omega$	±0.0001 Ω/°C	0.0001Ω	
	0 to 100Ω	2 mA	$\pm 0.008\Omega$	±0.0002 Ω/°C	0.0002Ω	
	0 to 300Ω	1 mA	$\pm 0.015\Omega$	±0.0004 Ω/°C	0.0003Ω	
	0 to $1 k\Omega$	500 μΑ	$\pm 0.05\Omega$	±0.001 Ω/°C	0.0007Ω	
	0 to $3 k\Omega$	200 µA	$\pm 0.1\Omega$	±0.003 Ω/°C	0.002Ω	
	0 to $10\mathrm{k}\Omega$	50 μA	$\pm 0.25\Omega$	±0.01 Ω/°C	0.007Ω	
	$0 \text{ to } 30 \text{ k}\Omega$	50 μA	±1 Ω	$\pm 0.02 \Omega/^{\circ}C$	0.008Ω	
	0 to $100\mathrm{k}\Omega$	5 μA	$\pm 4\Omega$	±1 Ω/°C	0.12Ω	
	0 to $300\mathrm{k}\Omega$	5 μΑ	$\pm 13 \Omega$	±2Ω/°C	0.2Ω	
	0 to $2.5\mathrm{M}\Omega$	1 μA	$\pm 1 \mathrm{k}\Omega$	±50 Ω/°C	10Ω	
Thermistors	0 to 10Ω	1 mA	$\pm 0.007\Omega$	$\pm 0.0002 \Omega/^{\circ}C$	0.0003Ω	
	0 to 30Ω	300 μΑ	$\pm 0.03\Omega$	$\pm 0.0004 \Omega/^{\circ}C$	0.001Ω	
	0 to 100Ω	100 μA	$\pm 0.07\Omega$	$\pm 0.002\Omega/^{\circ}C$	0.002Ω	
	0 to 300Ω	30 μΑ	$\pm 0.25\Omega$	$\pm 0.004 \Omega/^{\circ}C$	0.006Ω	
	0 to 1 k Ω	10 μΑ	$\pm 0.6\Omega$	±0.01 Ω/°C	0.02Ω	
	0 to $3 k\Omega$	3 μΑ	$\pm 2\Omega$	$\pm 0.06\Omega/^{\circ}\mathrm{C}$	0.06Ω	
	$0 \text{ to } 10 \text{ k}\Omega$	1 μA	$\pm 6\Omega$	$\pm 0.2\Omega/^{\circ}C$	0.2Ω	
	$0 \text{ to } 30 \text{ k}\Omega$	300 nA	$\pm 25\Omega$	±1 Ω/°C	1.0Ω	
	0 to $100 \mathrm{k}\Omega$	100 nA	$\pm 150\Omega$	±3 Ω/°C	6Ω	
	0 to $300\mathrm{k}\Omega$	30 nA	$\pm 1\mathrm{k}\Omega$	$\pm 20\Omega/^{\circ}C$	40Ω	
	0 to $2.5\mathrm{M}\Omega$	1 μΑ	$\pm 1\mathrm{k}\Omega$	±30 Ω/°C	10Ω	
		•				

Ordering Information

CTC100 Cryogenic temperature controller \$2600
Option 01 GPIB interface (replaces RS-232) \$495
O100CTRM Rack mount tray \$150



phone: (408)744-9040 www.thinkSRS.com