iSource+™ Rubidium Clock Lifetime and MTBF Data

Spectratime

iPrecision Timing Solutions®

November 2014

LPFRS / RMO / STARLPRO / LNRClok / GRClok Rubidium Clock Lifetime and MTBF

Contents	
Life Time	1
MTBF CalculationLPFRS / RMO	
StarLPRO / LNRClok / GRClok	3

Life Time

The aging of the Rb lamps due to the diffusion of rubidium metal into the glass walls of its bulbs, has been previously studied by Mr, C.-H. Volk and Mr, R.-P. Frueholz, of Aerospace Corporation, CA, USA.

These studies have been carried out using the differential scanning calorimetry technique.

The calorimetric technique measures the elemental Rubidium content in a lamp used in normal conditions.

TNT is using the SCHOTT glass 8436 type, which has been proved to be presently the best alkali resistant glass on the market .

The figure 1 shows the rubidium consumption of a lamp made of SCHOTT 8436 glass, and shows a logarithmic curve of the rubidium consumption, which is approximately 100 µg after a period of one year of operation.

Together with the Observatory of Neuchâtel, TNT performed the same measurement on one sample lamp n° 2285, in order to confirm the Aerospace Corporation data. A very similar characteristics was obtained (see figure 2).




Fig 1: Aerospace corp. data

These Rb consumption characteristics vary according to the square root of time:

$$M(t) = (0.9 \pm 0.15) \times \sqrt{t}$$

where: M=Rb mass in mg t=time, in hours

After 20 years of operation, $M = 376 \mu g$.

The rubidium lamps manufactured by T.N.T. are filled with 700 µg to 1000 µg of rubidium.




Fig 2: ON/TNT data

Consequently, the average life time, according to the diffusion of the rubidium into the glass walls, is:

$$t[h] = \left(\frac{700}{1.05}\right)^2 = 444'444 \text{ hours (50 years)}$$

MTBF Calculation

LPFRS / RMO

System MTBF Report									
System Data Records: 3			File Rev : 6	Model: Active K-of-N		MIL-HDBK-217F			
Miniature Rubidium clock Basic FR			Basic FR= 4.0066	066 f / million hours MTBF			= 249585.1260 hours		
Rec#	DC	Qty	TFR [%]	FR	File rev	Env	°C	Parts	
1	Circuit: CELLASSY			LPFRS/RMO CELL ASSEMBLY					
	100%	1	0.025050 [0.63%]	0.025050	5	GB	25°C	3	
2	Circuit: LAMPASSY			LPFRS/RMO LAMP ASSEMBLY					
	100%	1	0.027194 [0.68%]	0.027194	6	GB	25°C	8	
3	Circuit: MAIN BOARD			LPFRS/RMO MAIN BOARD					
	100%	1	3.9544 [98.7%]	3.9544	36	GB	25°C	152	

StarLPRO / LNRClok / GRClok

System MTBF Report								
System Data Records: 3			File Rev : 1	Model: serial		Bellcore	Bellcore 6	
Miniature Rubidium clock Basic			Basic FR= 2.5096	MTBF= 398'461.646 hours				
Rec#	DC	Qty	TFR [%]	FR	File rev	Env	°C	Parts
1	Circuit: MAIN BOARD			StarLPRO / LNRClok / GRClok MAIN BOARD				
	100%	1	2.169 [86.45%]	2.169	1	GB	25°C	333
2	Circuit: CELL and LAMP			StarLPRO / LNRClok / GRClok CELL and LAMP ASSEMBLY				
	100%	1	0.216 [8.61%]	0.216	1	GB	25°C	23
3	Circuit: OSCILLATOR			StarLPRO / LNRClok / GRClok OSCILLATOR ASSEMBLY				
	100%	1	0.124 [4.94%]	0.124	1	GB	25°C	39